Enzymax, LLC

Product Information

Product Name:EZ DNA Gel Extraction KitRelated products:Catalog No.:EZGE301DNA mini spin columnSize:50 prepsDNA Tini spin columns

Related products: DNA mini spin columns (30-50μl vol.) DNA Tini spin columns (5-20μl vol.) EZ DNA Clean-up kit (50)

Catalog No. EZC101 EZC106 EZDC 201

Features:

- Easy and rapid with 20 min procedure using spin column

- High DNA recovery

-Extract and purify DNA fragments (70 bp-12 Kb) from standard or low-melting agarose gel in TAE or TBE buffer - It can be used for easy and rapid DNA cleanup from mini prep., PCR amplification (from 80 bp-12 Kb), and enzymatic reactions (e.g., labeling, restriction, and dephosphorylation etc.)

-Purified DNA can be used directly in applications such as restriction, ligation, transformation, PCR amplification, labeling, hybridization, and sequencing

Kit Contents:

Components	EZ DNA Gel extraction 50 Preps (cat# EZGE301)
Gel Extraction Buffer 5xWash Buffer*	25 ml 5 ml
wini spin column with collection tube	50

*Add 100% ethanol before use: add 20 ml 100% ethanol to 5 ml 5xWash buffer.

Caution:

Gel Extracion Buffer contains chaotropic salt. Please use proper safety precautions and always wear gloves when handling the reagent. Avoid contact with skin, eyes or clothing. In case of accidental spill or contact, wash thoroughly with water, seek medical advice if necessary.

Protocols:

This kit is designed for DNA fragment recovery from agarose gel using Mini spin columns (Cat#EZC101). If you have small amount of DNA sample on the gel, you can use Tini Spin Columns (EZC106) with same protocol but cut down all the solution to $\frac{1}{2}$ to $\frac{1}{3}$. The elution volume for Tini spin column can be as low as 5μ l.

1. Add 3 volumes of the Gel Extracion Buffer into gel slice.

2. Incubate at 65°C for about 10 min till the gel slice is completely dissolved. Increase the temperature to 85°C, incubation time, or add more extraction buffer if the gel concentration is more than 2%. <u>Note:</u> If the color of the mixture turns a blue or purple color, adjust pH by adding a small volume of 3M Sodium acetate (pH 5.0).

3. Load the sample mixture onto the **Mini Spin Column** (or Tini spin Column) and spin in a microcentrifuge for 1 min at full speed (about 10,000 rpm). Do not load more than 700 μ l of sample on the Mini spin column (or 400 μ l on Tini spin column) at one time. Discard the flow through and load more sample mixture if needed.

4. Wash the column by adding 400 µl of 1x Wash Buffer (ethanol added) and centrifuge for 1 min.

5. Wash once by adding 500 µl of 80% ethanol and centrifuge for 1 min.

6. Discard flow through and place the column back in the same tube and spin for 1 min. You can also cut the cap off (this will help remove ethanol completely).

7. Place the column in a clean 1.5 ml micro-centrifuge tube.

8. Add 30 μ l or more **distilled water** or **TE** (10 mM Tris-HCl, 1 mM EDTA, pH 8.0) (preheated at 65°C for better yield) to the **center** of the column and leave at room temperature for 5 min. Spin the column for 1 min to elute DNA from the column. For Tini spin column, add as low as 5 μ l **distilled water** or **TE** for elution.

FOR RESEARCH USE ONLY

Enzymax LLC, 870 Corporate Dr. Suite 201, Lexington, KY 40503 Tel. (859) 219-8482; Fax. (859) 219-0653; Web: www.enzymax.net