


Enzymax, LLC

A345 AStECC Building, University of Kentucky, Lexington, KY 40506
Tel. 859-257-2300 ext. 280; Fax: 859-257-2489; www.enzymax.net

Unstained Protein Marker, Broad Range I (10-200kDa), 2x250ul for 100 mini gel applications, 5ul/well.

This unstained protein marker contains integral Strep-tag II sequence. This sequence can be detected directly on Western blots using a reporter enzyme conjugated with Strep-Tactin or Strep-tag II specific antibody.



8-16% SDS-PAGE gel Stained with Coomassie blue

8-16% SDS-PAGE gel blot detected with Strep-Tactin-AP conjugate