


Enzymax, LLC

A345 AStEC Building, University of Kentucky, Lexington, KY 40506
Tel. 859-257-2300 ext. 280; Fax: 859-257-2489; www.enzymax.net

1 KB plus DNA Marker (75bp-20Kb), quantitative and ready to use.



Range

15 fragments (in bp): 20000, 10000, 7000, **5000**, 4000, 3000, 2000, **1500**, 1000, 700, **500**, 400, 300, 200, 75.

Molecular Biology
FOR RESEARCH USE ONLY